

## Product Description

An Aerosol Mask is a device placed over a patient nose or mouth to administer aerosolized humidification or medication for patients to inhale.

Aerosol Therapy is a technique of administering medication directly into the airway and lungs

An aerosol is a suspension of liquid and/or solid particles.

Saudi Mais Co. for Medical Products

Third Industrial Area, Al-Kharj Road, Street No.256 P.O Box: 3900, Riyadh 14335-7599 Saudi Arabia



## Aerosol Mask





## Intended Use:

The aerosol mask is with a nebulizer or 22mm corrugated tubing for combined oxygen therapy and humidification.

The Primary use of aerosol mask therapy is treatment of respiratory disorders that include:

1. Obstructive lung diseases such as:
  - a. Asthma
2. Chronic Obstructive pulmonary disease
  - a. Bronchitis
  - b. Emphysema
3. Bronchiectasis

## Size:

Ref. Code	Description
341-1000	Pediatric / Small
341-2000	Child / Medium
341-3000	Adult / Large

## Technical Features

1. Soft and Comfortable face mask.
2. Adjustable elastic Strap.
3. Standard 22mm Male Connector.
4. Can be used with nebulizer or 22mm corrugated tube.
5. For combined oxygen therapy and humidification.
6. Made of medical grade material.
7. Non-Sterile

## Benefits and Features:

- Individually packaged.
- Direct delivery to the treatment site.
- Faster onset action than oral medication
- Lower requirement of dosage than systemic administration.

## Contraindications:

- Aerosol therapy is contraindicated when there is a known: Hypersensitivity or history of an allergic reaction to a specific pharmacological agent and / or its excipients

## Warnings & Precautions

1. For single patient use.
2. Please authenticate the integrity of all parts and pack-aging before use. If there is damage, please do not use this product.
3. Dispose of as specified by local medical waste disposal regulations.
4. Any modification, rework, or sterilization and will in-crease the risk of cross-infection of disease.
5. Discontinue the use if it causes discomfort and consult with physician.
6. Only used by health care professionals.

## Contact Us

Saudi Mais Co. for Medical Products  
Third industrial Area, Street No.: 256  
P.O Box:3900, Riyadh 14335 - 7599

+966 11 265 0184  
[medical@mais.com.sa](mailto:medical@mais.com.sa)

Visit us on the Web:  
[www.maisco.com.sa](http://www.maisco.com.sa)