


Product Description

A Foley Catheter is a flexible tube that a clinician passes through the urethra and into the bladder to drain urine.

It is the most common type of indwelling urinary catheter.

Indwelling urinary catheters are most commonly used to assist patients who cannot urinate on their own.



Foley Catheter

Saudi Mais Co. for Medical Products

Third Industrial Area, Al-Kharj Road, Street No.256 P.O Box: 3900, Riyadh 14335-7599 Saudi Arabia





Size & Specifications:

2-Way Foley Balloon Catheter				
Cat.Ref	Size CH	Balloon Capacity ml	Length cm	
442-0630	06	3-5	30	
442-0830	08	3-5	30	
442-1030	10	5-15	30	
442-1240	12	5-15	40	
442-1440	14	30-50	40	
442-1640	16	30-50	40	
442-1840	18	30-50	40	
442-2040	20	30-50	40	
442-2240	22	30-50	40	
442-2440	24	30-50	40	
442-2640	26	30-50	40	
3-Way Foley Balloon Catheter				
Cat.Ref	Size CH	Balloon Capacity ml	Length cm	
443-1440	14	30-50	40	
443-1640	16	30-50	40	
443-1840	18	30-50	40	
443-2040	20	30-50	40	
443-2240	22	30-50	40	
443-2440	24	30-50	40	
443-2640	26	30-50	40	
All Silicone 2-way Foley Balloon Catheter				
Cat.Ref 2-way	Cat.Ref 3-way	Size CH	Balloon Capacity ml	Length cm
444-0630	--	06	3-5	30
444-0830	--	08	3-5	30
444-1030	--	10	5-15	30
444-1240	--	12	5-15	40
W-1440	445-1440	14	30-50	40
444-1640	445-1640	16	30-50	40
444-1840	445-1840	18	30-50	40
444-2040	445-2040	20	30-50	40

Intended Use:

Foley Catheters are used in the following conditions:

1. On patients who are anesthetized or sedated for surgery or other medical care.
2. On comatose patients.
3. On some incontinent patients.
4. On patients whose prostate is enlarged to the point that urine flow the bladder is cut off.
5. On patients with acute urinary retention.

Technical Features:

- Made of 100% silicon coated natural latex.
- Symmetrical & safe balloon inflation.
- Soft & smooth rounded closed tip.
- Two lateral eyes.
- Color coded.


Complications:

There are several risks in using Foley catheter:

1. The Balloon might not inflate after it is in place.
2. Urine stops flowing into the bag.
3. The urethra begins to bleed

Warnings:

- Only use by professional healthcare practitioner.
- Do not exceed recommended capacities.
- Designed for Single use Only.
- Sterility will not be guaranteed if the package is damaged / opened.
- Do not resterilize any portion of the device.



Contact Us

Saudi Mais Co. for Medical Products
Third industrial Area, Street No.: 256
P.O Box:3900, Riyadh 14335 - 7599

+966 11 265 0184

medical@mais.com.sa

Visit us on the Web:

www.maisco.com.sa